

endo★star **PEESO REAMERS****1. Description:**

Endostar Peeso Reamers drills are engine-driven instruments and they are used to enlarge root canal orifice and coronal third of the canals. Peeso reamers differ from Gates Glidden drills in that the blades spread over a wide surface and the shape is cylindrical.

- slightly flexible,
- rotation speed: 800-1200 rpm,
- non cutting tip.

2. Recommended number of use:

Maximum of 5-10 times, depending on the size (see Table 1), provided that visual inspection performed by the practitioner prior to use shows that the instrument remains undamaged, is not bent, deformed, does not show signs of blade wear. Prolonging the life of the instrument more than recommended may result in the blade breaking.

Table 1

File No	01	02	03	04	05	06
Number of times that the instrument can be used	5	5	5	10	10	10

3. Cleaning and disinfection:

Detailed instructions for cleaning and disinfection can be found at www.poldent.pl and www.endostar.eu.

4. Sterilization:








This is a non-sterile product. Sterilize before use. The instruments can be sterilized in a steam sterilizer (autoclave) at 134°C. Recommended sterilization time 35-40 minutes. Instruments can be disinfected with mild disinfectants and washed in ultrasonic cleaners.

5. Storage:

Instruments should be stored at room temperature in a dry, dust-free and clean environment.

6. Warnings:

Only for dental use.

	CE mark and identification number of the notified body
	Sterilize in a steam autoclave at 134 ° C
	Non-sterile product
	Used for root canal preparation.
	Serial number/ Lot number
	Consult instruction for use
	Stainless steel



Manufactured by:

Poldent Co. Ltd.

Dzika 2 Street

00-194 Warsaw, POLAND

tel.: +48 22 351-76-50 do 51

fax: +48 22 351-76-79

www.poldent.pl, www.endostar.eu, e-mail: poldent@poldent.pl, endostar@endostar.eu



This manual was last updated: 01.2018