

## PMMA Premia:

### PMMA Forms:

'98 mm (Wieland-comp.)

'95 mm (ZirkonZahn-comp.)

'71 mm (Ceramill-comp.)

### Shades:

A1 A2 A3 B1 B2 B3;

Bleach: BL1, BL2, BL3;

Clear.

12, 14, 15, 16, 18, 20, 25 mm

Multilayer PMMA: A1, A2, A3, B1.

15, 18, 20, 25 mm

- ▶ outstanding millability
- ▶ precise colour match
- ▶ bending strength: 114 MPa
- ▶ Vickers hardness: 26.6 HV
- ▶ E-Module: 2771 MPa
- ▶ residual monomer < 1%

## Features:

- ✓ Colours & translucent.
- ✓ Optimal surface quality, reproducing fissures. Antagonist-friendly.
- ✓ Temporary restorations fabricated easily.
- ✓ Venerable with acrylic veneering resins.
- ✓ Biocompatible: void of irritant substances.
- ✓ To avoid melting during milling & and bur filled up with residue, use low spindle RPM and low federate!

## Specs:

Time in mouth up to 5 yrs.

Full contour: ✓

Framework: ✓

Gingival cover: ✓

Colour accuracy: **excellent.**

Millability: **excellent.**

Flex. strength: 114 MPa.

Water: insoluble.

Residual monomer: < 1,0%. (100% PMMA).